



# Recombinant Human Upstream-binding protein 1 (UBP1)

<b>Product Code</b>	CSB-EP865191HU
<b>Abbreviation</b>	UBP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NZI7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAWVLKMDEV IESGLVHDFD ASLSGIGQEL GAGAYSMSDV LALPIFKQED SSLPLDGETE HPPFQYVMCA ATSPAVKLHD ETLTYLNQGQ SYEIRMLDNR KMGDMPEING KLVKSIIRVV FHDRRLQYTE HQQLEGWKWN RPDRLDLD IPMSVGIIDT RTNPSQLNAV EFLWDPKRT SAFIQVHCIS TEFTPRKHGG EKGVPFRIQV DTFKQENGE YTDHLHSASC QIKVFKPKGA DRKQKTDREK MEKRTAHEKE KYQPSYDTTI LTEMRLPII EDAVEHEQKK SSKRTLPAKY GDSLAKRGSC SPWPDAPTAY VNNSPSPAPT FTSPQQSTCS VPDSNSSSPN HQGDGASQTS GEQIQPSATI QETQQWLLKN RFSSYTRLFS NFSGADLLKL TKEDLVQICG AADGIRLYNS LKSRSVRPRL TIYVCREQPS STVLQGGQQA ASSASENGSG APYVYHAIYL EEMIASEVAR KLALVFNIPL HQINQVYRQG PTGIHILVSD QMVQNFQDES CFLFSTVKAE SSDGIHILK
<b>Source</b>	E.coli
<b>Target Names</b>	UBP1
<b>Protein Names</b>	Recommended name: Upstream-binding protein 1 Alternative name(s): Transcription factor LBP-1
<b>Expression Region</b>	1-540
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.