



# Recombinant Human CCR4-NOT transcription complex subunit 2 (CNOT2)

<b>Product Code</b>	CSB-EP865193HU-B
<b>Abbreviation</b>	CNOT2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NZN8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVRTDGHTLS EKRNYQVTNS MFGASRKKFV EGVDSYHDE NMYYSQSSMF PHRSEKDMLA SPSTSGQLSQ FGASLYGQQS ALGLPMRGMS NNTPQLNRSL SQGTQLPSHV TPTTGVPTMS LHTPPSPSRG ILPMNPRNMM NHSQVGQGIG IPSRTNSMSS SGLGSPNRSS PSIICMPKQQ PSRQPFTVNS MSGFGMNRNQ AFGMNNLSLSS NIFNGTDGSE NVTGLDLSDF PALADRNRRE GSGNPTPLIN PLAGRAPYVG MVTKPANEQS QDFSIHNEDF PALPGSSYKD PTSSNDDSKS NLNTSGKTTS STDGPKFPGD KSSTTQNNNQ QKKGIVLPD GRVTNIPQGM VTDQFGMIGL LTFIRAAETD PGMVHLALGS DLTTLGLNLN SPENLYPKFA SPWASSPCRP QDIDFHVPSE YLTNIHIRDK LAAIKLGRYG EDLLFYLYYM NGGDVLQLLA AVELFNRDWR YHKEERVWIT RAPGMEPTMK TNTYERGTYY FFDCLNWRKV AKEFHLEYDK LEERPHLPST FNYNPAQQAF
<b>Source</b>	E.coli
<b>Target Names</b>	CNOT2
<b>Protein Names</b>	Recommended name: CCR4-NOT transcription complex subunit 2 Alternative name(s): CCR4-associated factor 2
<b>Expression Region</b>	1-540
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.