



Recombinant Human DnaJ homolog subfamily C member 27 (DNAJC27)

Product Code	CSB-BP865194HU
Abbreviation	DNAJC27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZQ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEANMPKRKE PGRSLRIKVI SMGNAEVGKS CIIKRYCEKR FVSKYLATIG IDYGVTKVHV RDREIKVNIF DMAGHPFFYE VRNEFYKDTQ GVILVYDVGQ KDSFDALDAW LAEMKQELGP HGNMENIIFV VCANKIDCTK HRCVDESEGR LWAESKGFLY FETSAQTGEG INEMFQTFYI SIVDLCENGG KRPTTNSSAS FTKEQADAIR RIRNSKDSWD MLGVKPGASR DEVNKAYRKL AVLLHPDKCV APGSEDAFKA VVNARTALLK NIK
Source	Baculovirus
Target Names	DNAJC27
Protein Names	Recommended name: DnaJ homolog subfamily C member 27 Alternative name(s): Rab and DnaJ domain-containing protein
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.