



Recombinant Human Homologous-pairing protein 2 homolog (PSMC3IP)

Product Code	CSB-MP865213HU
Abbreviation	PSMC3IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P2W1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSKGRAEAAA GAAGILLRYL QEQRNPYSSQ DVFGNLQREH GLGKAVVVKT LEQLAQQGKI KEKMYGKQKI YFADQDQFDM VSDADLQVLD GKIVAL TAKV QSLQQSCRYM EAEKELSSA LTTPEMQKEI QELKKECAGY RERLKNIAA TNHVTPEEKE QVYRERQKYC KEWRKRKRMA TELSDAILEG YPKSKKQFFE EVGIETDEDY NVTLPDP
Source	Mammalian cell
Target Names	PSMC3IP
Protein Names	Recommended name: Homologous-pairing protein 2 homolog Alternative name(s): Nuclear receptor coactivator GT198 PSMC3-interacting protein Proteasome 26S ATPase subunit 3-interacting protein Tat-binding protein 1-interacting protein
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.