



Recombinant Schizosaccharomyces pombe Transcription factor IIIB 60 kDa subunit (brf1)

Product Code	CSB-BP865244SXV
Abbreviation	brf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P6R0
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MGCPNCGSTT FESDTASGNT YCTQCGVVVE QDAIVSEVTF GEASTGAAVV QGSLVSNDQT HARTFGGPYR NQGSVESREL TIANGRRRIS ALAIALKLNE RHIEAAVRYF TLAINNNFIK GRRSQYVVAS CLYIVCRISK TSHMLIDFSD ILQINVFKLG STFLKLCRVL RPNLPLLDPS LYISRFASLL EFGPETHRVA NDAIRLVARM NRDWMQIGRR PAGICGACLL IAARMNNFRR SVREVVHVVK VADITIQKRL DEFKLTESGD LSIADFRNIW LEGQSDPPSF TKNQKFQQYG AQKVSNIIDHT QEYMSPIKRT PDFDGENVKS EELSQTVKVE SQETPVHLKA DEREIRKEVT ETLKGDELK ISLQVNVKFS EEEVTLEDVD DDEIEDILLD KDEILTKTQV WMELNKDYLA EEEAKNLKLQ EDLKKGIVRQ PRKRRRRYRPR DSTSDGIADT AAESAKEMMQ QRAFSKKINY EALDMLFDEE QS
Source	Baculovirus
Target Names	brf1
Protein Names	Recommended name: Transcription factor IIIB 60 kDa subunit Short name= TFIIIB Alternative name(s): B-related factor 1 Short name= BRF-1 TFIIIB-related factor
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.