



# Recombinant Schizosaccharomyces pombe Transcription factor IIIB 60 kDa subunit (brf1)

<b>Product Code</b>	CSB-EP865244SXV-B
<b>Abbreviation</b>	brf1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P6R0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGCPNCGSTT FESDTASGNT YCTQCGVVVE QDAIVSEVTF GEASTGAAVV QGSLVSNQDT HARTFGGPYR NQGSVESREL TIANGRRRIS ALAIALKLNE RHIEAAVRYF TLAINNNFIK GRRSQYVVAS CLYIVCRISK TSHMLIDFSD ILQINVFKLG STFLKLCRVL RPNLPLLDPS LYISRFASLL EFGPETHRVA NDAIRLVARM NRDWMQIGRR PAGICGACLL IAARMNNFRR SVREVVHVVK VADITIQKRL DEFKLTESGD LSIADFRNIW LEGQSDPPSF TKNQKFQQYG AQKVSNIHT QEYMSPIKRT PDFDGENVKS EELSQTVKVE SQETPVHLKA DEREIRKEVT ETLKGDELK ISLQVNVKFS EEEVTLEDVD DDEIEDILLD KDEILTKTQV WMELNKDYLA EEEAKNLKLQ EDLKKGIVRQ PRKRRRRYRPR DSTSDGIADT AAESAKEMMQ QRAFSKKINY EALDMLFDEE QS
<b>Source</b>	E.coli
<b>Target Names</b>	brf1
<b>Protein Names</b>	Recommended name: Transcription factor IIIB 60 kDa subunit Short name= TFIIIB Alternative name(s): B-related factor 1 Short name= BRF-1 TFIIIB-related factor
<b>Expression Region</b>	1-492
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.