



Recombinant Schizosaccharomyces pombe ZZ-type zinc finger-containing protein P35G2.11c (SPBP35G2.11c)

Product Code	CSB-BP865250SXV
Abbreviation	SPBP35G2.11c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P792
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MLCNRTCPSS DCGILCNHSV PSFPPFHSSV ANIHFTKENN LKSNIFEHNN NSPTLRSSSV ACNTCLKIIR NDSFHCTKCF DFDVCRDCYA KQAFHLPCPK PHFVLVRSSI PSVASLTCSV NSMSVSPQSN FMYAICDHCE QPIHNVRYKC SVCDDYDICE SCLTDNSHSN THAFVRITKA YPHGLPSFHL FPQFLSAELF ESASTVHRSV QCDNCLAHPI VGPRFHCLVC EDYDLCSSCV SHVHHDHHS LRLTREISAS PLHLSKPEKP KVLNDFDKLV EDSILPLELS PGCPFYKIWH IRNTSCQSWP SPLYVKFNGG DKLFPGDNPY SFPITSSVHP GEDVNFTVAL KVPEKSNKEI FTAFFNICSD DGSVFHKS LC AFLRVPKSFS RISKSSF
Source	Baculovirus
Target Names	SPBP35G2.11c
Protein Names	Recommended name: ZZ-type zinc finger-containing protein P35G2.11c
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.