



Recombinant Danio rerio Homeobox protein Hox-B10a (hoxb10a)

Product Code	CSB-MP865529DIL
Abbreviation	hoxb10a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9PWL7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MALHRHPFVH GAASWDGEQP SPVQLPHISA CPFTSSGRKE DPFYFTLDPT GQARQPLDIS AFSRIMTDMG SLNSADDHRP ETSFYSGHKL LSLNTDTPDPE SPSLSSHSDS QHLHSLLSLSA PCSETDNKHD MQYLAMESTK YPSQWSESRT NRSIITLSI SQRSKDDIEP>NNLQTDFTRG DKTPREKTQD VTLENAANGW LSAKAGRKKR CPYSKHQILE LEKEFLFNMV LTRERRLEIS RSINLTDRQV KIWFQNRMMK LKKMTREHRT RDPGTSFTV
Source	Mammalian cell
Target Names	hoxb10a
Protein Names	Recommended name: Homeobox protein Hox-B10a Short name= Hox-B10
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.