



Recombinant Tomato mosaic virus Capsid protein (CP)

Product Code	CSB-EP865535TGG
Abbreviation	CP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Q1T6
Product Type	Recombinant Protein
Immunogen Species	Tomato mosaic virus (strain Kazakh K1) (ToMV) (TMV strain K1)
Purity	>85% (SDS-PAGE)
Sequence	SYPITSPSQ FVFLSSVWAD PIELLNICTN SLGNQFQTQQ ARTTVQQQFS EWWPEFPQST VRFPGDVYKV YRYNAVLDPD ITALLGTFDT RNRRIEVENR QSPTTAETLD ATRRVDDATV AIRSAINNLV NELVRGTGLY NQNTFESMSG LVWTSAPAS
Source	E.coli
Target Names	CP
Protein Names	Recommended name: Capsid protein Alternative name(s): Coat protein
Expression Region	2-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.