



Recombinant Mouse Dickkopf-related protein 3 (Dkk3)

Product Code	CSB-MP865568MO
Abbreviation	Dkk3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9QUN9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	PSPTVTWT PAEPGPALNY PQEEATLNEM FREVEELMED TQHKLRSAVE EMEAEAAAK TSSEVNLASL PPNYHNETST ETRVGNNTVH VHQEVHKITN NQSGQVVFSE TVITSVGDEE GKRSHEDIID EDCGPTRYCQ FSSFKYTCQP CRDQQMLCTR DSECCGDQLC AWGHCTQKAT KGGNGTICDN QRDCQPGLCC AFQRGLLFPV CTPLPVEGEL CHDPTSQLLD LITWELEPEG ALDRPCASG LLCQPHSHSL VYMCKPAFVG SHDHSEESQL PREAPDEYED VGFIGEVQRQ LEDLERSLAQ EMAFEGPAPV ESLGEEEE
Source	Mammalian cell
Target Names	Dkk3
Protein Names	Recommended name: Dickkopf-related protein 3 Short name= Dickkopf-3 Short name= Dkk-3 Short name= mDkk-3
Expression Region	23-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. The expression of this gene is decreased in a variety of cancer cell lines and it may function as a tumor suppressor gene. Alternative splicing results in multiple transcript variants encoding the same protein.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.