



Recombinant Mouse LIM domain-containing protein 1 (Limd1)

Product Code	CSB-EP865576MO
Abbreviation	Limd1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9QXD8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MDKYDDLGLE ASKFIEDLNM YEASKDGLFR VDKGAGNPE FEETRRVFAT KMAKIHLLQQQ QQQQLLQEEA LPRAGRSPVN GGNRQGASGK LAADGAAKPP LAVPTVAPGL ATTTAAQPS YPSQEQRIRP SAHGARPGSQ NCGSREGPVS SQRPALHGLS PSCEDPSCLT HGDYYDNFSL ASPQWGDKPE GCPSVSLGVG SGWPGCPGND STLPKSCGDH HPYQPQLSTV CSGRSFESGI SGQDGGIGGH SSEKPTGLWS TASSQRVNLG FSSMGLENGT SAQPKGTTVS APMVPSSASQ GACPKRDSGL GYEASGRVFK PLVDTQPWLQ DGPKSYLSVS APLSSTAGKD STQPGMTTGL DPKFGCVESG TSPKPSPTSN VHPVMSTPSE LSCKESSPSW STDSSLEPVL PGSPTPSRVR LPCQTLAPGP ELGPSTAEK LEALTQRLER EMDAHPKADY FGSCVKCSKG VFGAGQACQA MGDLYHDACF TCAACSRKLR GKAFYFVNGK VFCEEDFLYS GFQQSADRCF LCGHLIMDMI LQALGKSYHP GCFRCVICNE CLDGVPFTVD SENKIYCVRD YHKVLAPKCA ACGLPILPPE GSDETIRVVS MDRDYHVECY HCEDCGLELN DEDGHCRCYPL EDHLFCHSCH VKRLEKGPSP APLHQHHF
Source	E.coli
Target Names	Limd1
Protein Names	Recommended name: LIM domain-containing protein 1
Expression Region	1-668
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.