



# Recombinant Mouse Protein DBF4 homolog A (Dbf4)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP865604MO-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Abbreviation</b>      | Dbf4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | Q9QZ41                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | MNLETMRIHS KAPLPGGIQD RNEKNRPSLK SLKADNRLEK SKYKPLWGKI<br>FYLDLPSITI CEKLQKDIKE LGGRVEEFLS KDISYFVSNK KEAKYAQTLG<br>RVSPVPSPES AYTAETTSPH PSHDGSSFKS QDRVCLSRGK LLAEKAVKDH<br>DFIPANSILS NALSWG VKIL HIDDIRYYIE QKKKALSALK KSSASGKDAG<br>KKAGPGIQKT RTGRLKKPFL KVEDVNRCYR PFYLQLPSLP CINYFLQKPC<br>SPFDIEKSSS VQKQAQPKLR INMDGDKCGT PVQLQLKEKR KKG YCECCLQ<br>KYEDLETHLL SEKHRNFAQS NQYQVDDIV SQLVDFVEY GRDTPQKKRI<br>RYSVGSLSV SANVLKNTAP KEKPLLEPNF QKDVGESESGH LLKPNSQYEE<br>TQKPEEKHGF ASEPTTYSSA GLKGCDRKPV SMFNASEPDP EQEYAQLPLR<br>DSTPEHQVTE GRNDGEQRVD PAPGVSQSCG QVSHLSTESN LPQPQLAADI<br>TQLSAKDLQE KGFHV VIGHA SDLVALNTSK EQLTMKARTP PCSPQEPHEC<br>DTENMENLPC GKIQRKVRML LGQQKANAEP SAELDKK RTE YLPAHEDRTC<br>GSPVQSLLDL FQTSEEKSEF LGFTGYTENS GICDVLDIWE EENSSTLLST<br>FFSSPSTSAF VGF |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Target Names</b>      | Dbf4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: Protein DBF4 homolog A Alternative name(s): MuDBF4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-663                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                                                                                                                        |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.