



Recombinant Mouse Protein DBF4 homolog A (Dbf4)

Product Code	CSB-MP865604MO
Abbreviation	Dbf4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9QZ41
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MNLETMRIHS KAPLPGGIQD RNEKNRPSLK SLKADNRLEK SKYKPLWGKI FYLDLPSITI CEKLQKDIKE LGGRVEEFLS KDISYFVSNK KEAKYAQTLG RVSPVPSPES AYTAETTSPPH PSHDGSSFKS QDRVCLSRGK LLAEKAVKDH DFIPANSILS NALSWGKIL HIDDRIYYIE QKKKALSALK KSSASGKDAG KKAGPGIQKT RTGRLKKPFL KVEDVNRCYR PFYLQLPSLP CINYFLQKPC SPFDIEKSSS VQKQAQPKLR INMDGDKCGT PVQLQLKEKR KKGyceCCLQ KYEDLETHLL SEKHRNFAQS NQYQVDDIV SQLVDFVEY GRDTPQKKRI RYSVGSLSV SANVLKNTAP KEKPLLEPNF QKDVGESESGH LLKPNSQYEE TQKPEEKHGF ASEPTTYSSA GLKGCDRKPV SMFNASEPDP EQEYAQLPLR DSTPEHQVTE GRNDGEQRVD PAPGVSQSCG QVSHLSTESN LPQPQLAADL TQLSAKDLQE KGFHVVGHA SDLVALNTSK EQLTMKARTP PCSPQEPHEC DTENMENLPC GKIQRKVRML LGQQKANAEP SAELDKKRT EYLP AHEDRTC GSPVQSLLDL FQTSEEKSEF LGFTGYTENS GICDVLDIWE EENSSTLLST FFSSPSTSAF VGF
Source	Mammalian cell
Target Names	Dbf4
Protein Names	Recommended name: Protein DBF4 homolog A Alternative name(s): MuDBF4
Expression Region	1-663
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.