



# Recombinant Mouse Protein DBF4 homolog A (Dbf4)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP865604MO                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Abbreviation</b>      | Dbf4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Uniprot No.</b>       | Q9QZ41                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | MNLETMRIHS KAPLPGGIQD RNEKNRPSLK SLKADNRLEK SKYKPLWGKI<br>FYLDLPSITI CEKLQKDIKE LGGRVEEFLS KDISYFVSNK KEAKYAQTLG<br>RVSPVPSPES AYTAETTSPPH PSHDGSSFKS QDRVCLSRGK LLAEKAVKDH<br>DFIPANSILS NALSWGKIL HIDDRIYYIE QKKKALSALK KSSASGKDAG<br>KKAGPGIQKT RTGRLKKPFL KVEDVNRCYR PFYLQLPSLP CINYFLQKPC<br>SPFDIEKSSS VQKQAQPKLR INMDGDKCGT PVQLQLKEKR KKGyceCCLQ<br>KYEDLETHLL SEKHRNFAQS NQYQVDDIV SQLVDFVEY GRDTPQKKRI<br>RYSVGSLSV SANVLKNTAP KEKPLLEPNF QKDVGESESGH LLKPNSQYEE<br>TQKPEEKHGF ASEPTTYSSA GLKGCDRKPV SMFNASEPDP EQEYAQLPLR<br>DSTPEHQVTE GRNDGEQRVD PAPGVSQSCG QVSHLSTESN LPQPQLAADI<br>TQLSAKDLQE KGFHVVGHA SDLVALNTSK EQLTMKARTP PCSPQEPHEC<br>DTENMENLPC GKIQRKVRML LGQQKANAEP SAELDKKRTY YLPAHEDRTC<br>GSPVQSLLDL FQTSEEKSEF LGFTGYTENS GICDVLDIWE EENSSTLLST<br>FFSSPSTSAF VGF |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Names</b>      | Dbf4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Names</b>     | Recommended name: Protein DBF4 homolog A Alternative name(s): MuDBF4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 1-663                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                                                                                                                   |



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.