



# Recombinant Histidine biosynthesis bifunctional protein HisB (hisB)

<b>Product Code</b>	CSB-MP865810EOD
<b>Abbreviation</b>	hisB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9S5G5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSQKYLFDLR DGTLISEPPS DFQVDRFDKL AFEPGVIPQL LKLQKAGYKL VMITNQDGLG TQSFPQADFD GPHNLMQIF TSQGVQFDEV LICPHLPADE CDCRKPVKL VERYLAEQAM DRANSYVIGD RATDIQLAEN MGINGLRYDR ETLNWPMIGE QLTRRDYAH VVRNTKETQI DVQVWLDREG GSKINTGVGF FDHMLDQIAT HGGFRMEINV KGDLYIDDHH TVEDTGLALG EALKIALGDK RGICRFGFVL PMDECLARCA LDISGRPHLE YKAEFTYQRV GDLSTEMIEH FFRSLSYTMG VTLHLKTKGK NDHHRVESLF KAFGRTLRQA IRVEGDTLPS SKGVL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	hisB
<b>Protein Names</b>	Recommended name: Histidine biosynthesis bifunctional protein HisB Including the following 2 domains: Histidinol-phosphatase EC= 3.1.3.15 Imidazoleglycerol-phosphate dehydratase Short name= IGPD EC= 4.2.1.19
<b>Expression Region</b>	1-355
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.