



Recombinant Arabidopsis thaliana Heat stress transcription factor A-8 (HSFA8)

Product Code	CSB-EP865820DOA-B
Abbreviation	HSFA8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9S7U5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVKSTDGGGG SSSSSSVAPF LRKCYDMVDD STTDSIISWS PSADNSFVIL DTTVFSVQLL PKYFKHSNFS SFIRQLNIYG FRKVDADRWE FANDGFVRGQ KDLLKNVIRR KNVQSSEQSK HESTSTTYAQ EKSGLWKEVD ILKGDQVLA QELIKVRQYQ EVTDTKMLHL EDRVQGMEES QQEMLSFLVM VMKNPSSLVQ LLQPKEKNTW RKAGEGAKIV EEVTDEGESN SYGLPLVTYQ PPSDNGTAK SNSNDVNDFL RNADMLKFCL DENHVPLIIP DLYDDGAWEK LLLLSPSRKK TKKQENIVKK GKDDLLEEE EEDGTMELDK SYMLKLISEE MEKPDDFEFG QLTPERSRNL EILTEQMELL ASNE
Source	E.coli
Target Names	HSFA8
Protein Names	Recommended name: Heat stress transcription factor A-8 Short name= AtHsfA8 Alternative name(s): AtHsf-03 Heat shock factor protein 5 Short name= HSF 5 Heat shock transcription factor 5 Short name= HSTF 5
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.