



Recombinant Arabidopsis thaliana F-box protein At3g49450 (At3g49450)

Product Code	CSB-MP865847DOA
Abbreviation	At3g49450
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SCL2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MMKKRKRKVS KENVALTISS SVGEYGENSG TLPTDLMVEI LSRVPAKSAA RFRCVSNLWN SLLRSPYLTN LFLKRSSARP HLLITFQAEG KWSFFSSPEY LISDQNSNLV VVDNHMDVVK DYSFGVCVPV CGLMCTSNEW VLSRKRDRM MICNPSTGQF KSLPKVRSCR GNVITYIGYN PIEKQYKVLK MTIREKPKFK KAEHQVLTG GTGKLRWML ECSVDHYPY HGSICINGVL FYLAMKSESK EYMTVSFHMK DENFMFIPNQ DLLSTLINYK GRLGGIRHKS FGFMDGGDVG FELWILDVVN QEWIRSIHVL PPMWKDVGGE TRVYFVGIIG SCEVVFSPFV KSNPFYIFHL DMKSNSITRV EIKGTGPLEG QAVYTFVNI ENVNLIM
Source	Mammalian cell
Target Names	At3g49450
Protein Names	Recommended name: F-box protein At3g49450
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.