



Recombinant Arabidopsis thaliana Acidic leucine-rich nuclear phosphoprotein 32-related protein (At3g50690)

Product Code	CSB-MP865849DOA
Abbreviation	At3g50690
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SCQ7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDEIWERAVE AALDGQTDRL ATRTLTLDGA VKCVQGR LPP PSVLEKFQNL QHLSVANIGV SSLEQFRLG NLQKLILSDN RITVGLLEFLV EAGLDSFCDL DLSNNRIQFV EDLAPLAELK LVSLDLYECP VTRLKDYRSR VFGLIKTLKY LDKTDAEGNE RPESDDEDDE EDEEDEEEEE EGDEEDPGSG EIDGDERAEA PRMSNGHSER VDGVDVDED EESDAEDDES EQATGVNGTS YRANGFRLEA VNGEEVREDD GDDSESGEEE VGEDNDVVEV HEIEDSENEE DGVDDEEDDE EDEEEEEVDN ADRGLGGSGS TSRLMNAGEI DGHEQGDDDE DGDGETGEDD QGVEDDGEFA DEDDDVEEED EESGEGYLVQ PVSQVEDHDA VGNDIEPINE DNDPDEEEEV EDDLPIPDQS LASSSRPKRK RDDDDDGEDD DDDDEER
Source	Mammalian cell
Target Names	At3g50690
Protein Names	Recommended name: Acidic leucine-rich nuclear phosphoprotein 32-related protein Alternative name(s): ANP32/acidic nuclear phosphoprotein-like protein
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.