



Recombinant Arabidopsis thaliana UDP-glycosyltransferase 83A1 (UGT83A1)

Product Code	CSB-EP865874DOA-B
Abbreviation	UGT83A1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SGA8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDNNSNKRMG RPHVVVIPYP AQQHVLPLIS FSRYLAKQGI QITFINTEFN HNRIISSLPN SPHEDYVGDQ INLVSIPDGL EDSPEERNIP GKLSESVLRF MPKKVEELIE RMMAETSGGT IICSVVADQS LGWAIEVAAK FGIRRTAFCP AAAASMVLGF SIQKLIDDGL IDSDGTVRVN KTIQLSPGMP KMETDKFVWV CLKNKESQKN IFQLMLQNNN SIESTDWLLC NSVHELETAA FGLGPNIVPI GPIGWASLE EGSTSLGSFL PHDRDCLDWL DRQIPGSVIY VAFGSFGVMG NPQLEELAIG LELTKRPVLW VTGDQQPIKL GSDRVKVVWR APQREVLSSG AIGCFVSHCG WNSTLEGAQN GIPFLCIPYF ADQFINKAYI CDVWKIGLGL ERDARGVVPR LEVKKKIDEI MRDGGEYEER AMKVKEIVMK SVAKDGISCE NLNKFVNWIK SQVN
Source	E.coli
Target Names	UGT83A1
Protein Names	Recommended name: UDP-glycosyltransferase 83A1 EC= 2.4.1.-
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.