



Recombinant Arabidopsis thaliana Putative Ras-related protein RABA4e (RABA4E)

Product Code	CSB-MP865910DOA
Abbreviation	RABA4E
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJZ5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAEEEVDYVF KVVINGDSAV GKSQLRARFT RDEFSMDSKA TIRCRFQYSN APYYKGAVGA MLVYDMTIRE SFEHIPQWLE ELRVHADKNI VIILIGNKTD LENQRSVPVE DAKEFAEKEG LFFLETSALN STNVENSFNT LLTEIFNKVN KKNLAKTTVS CSSQVSLLRP PCVAAQ
Source	Mammalian cell
Target Names	RABA4E
Protein Names	Recommended name: Putative Ras-related protein RABA4e Short name= AtRABA4e
Expression Region	1-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.