



Recombinant Arabidopsis thaliana Putative protein phosphatase 2C 53 (At4g08260)

Product Code	CSB-BP865984DOA
Abbreviation	At4g08260
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SUF4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEDRFSAITN LHGDHKQAIF GVVYVGHGGVK AAEFAAKNLD KNIVEEVVDA TFLKEEGFKG GSSCVTALVS EGSLVVSNA GDCRAVMSVGE MMNGKELKPR EDMLIRFTLW RIQGSLVVR GIGDAQLKKW VIAEPETKIS RVEHDHEFLI LASHGLWDKV SNQEAVDIAR PFCLRTEKPL LLAACKKLVD LSASRGSFDD ISVMLIPLRQ FV
Source	Baculovirus
Target Names	At4g08260
Protein Names	Recommended name: Putative protein phosphatase 2C 53 Short name= AtPP2C53 EC= 3.1.3.16 Alternative name(s): Protein phosphatase AP2C5
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.