



# Recombinant Arabidopsis thaliana Endoglucanase 18 (At4g09740)

<b>Product Code</b>	CSB-EP866029DOA
<b>Abbreviation</b>	At4g09740
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SZ90
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LDYGDALNK SILFFEGQRS GKLPTNQRVK WRADSGLSDG ASANVNLIGG YYDAGDNVKF VWPMSFTTTL LSWAALEYQN EITFVNQLGY LRSTIKWGTN FILRAHTSTN MLYTQVGDGN SDHSCWERPE DMDTPRTLYS ISSSSPGSEA AGEAAAALAA ASLVFKLVDS TYSSKLLNNA KSLFEFADKY RGSYQASCPF YCSHSGYQDE LLWAAAWLYK ATGEKSYLNY VISNKDWSKA INEFSWDNKF AGVQALLASE FYNGANDLEK FKTDVESFVC ALMPGSSSSQQ IKPTPGGILF IRDSSNLQYV TTATTILFY SKTLTKAGVG SIQCGSTQFT VSQIRNFAKS QVDYILGNNP LKMSYMGVGF TKYPTQPHHR GSSLPSIQSK PEKIDCNGGF SYYNFDTPNP NVHTGAIVGG PNSSDQYSDK RTDYSHAAPT TYINAAFIGS VAALISSS
<b>Source</b>	E.coli
<b>Target Names</b>	At4g09740
<b>Protein Names</b>	Recommended name: Endoglucanase 18 EC= 3.2.1.4 Alternative name(s): Endo-1,4-beta glucanase 18
<b>Expression Region</b>	22-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.