



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At4g11750 (At4g11750, At4g11745)

Product Code	CSB-MP866042DOA
Abbreviation	At4g11750, At4g11745
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9T0E2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MMNPPVDLPY IPDLLLLNCL ARVSRLYYPI LSLVSKRFRS LVASLELYEI RKLLGHREKC LYLNLRFSSSE SEPRWFTLCR RPTRSPNPNF NSGWFASCFI PHTMKKEKKS SDNNLMVPVP TSNFITPPSE LTRITNGSNI YIVRVFTNGA FSSRFLFMDC RSHTLHEVPR MDKTKKKPFM IRVLDGKIYV IEGCKNPDYS NLIERFDLKT QTWEHVPSPS AEIRGSYITG SLVYDGKLYL FGDKRVVYKP KENKWDVVGL EMPLRWTPSY ISCVVDNVIY SYGRSRVLMW YNTEERLWRY LKGLKLPKL PKDCTCVRLG YSGKVVVLW EKNVGGGDPQ KKMIWCAEIA LERRSANKIY GKIEWCDVVL TLPRSCSLLN FAAVTV
Source	Mammalian cell
Target Names	At4g11750
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At4g11750
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.