



Recombinant Uncharacterized protein Y116A8C.11 (Y116A8C.11)

Product Code	CSB-MP866171CXY
Abbreviation	Y116A8C.11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9U2T3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MESKQPEDIH QDCELKQDEL ISKLQHSQQN FDNIIKILAN FEKENENLAD ENEILKTRIA SLEKDVKQEK NATKMAKSLN EETINELTLS EDLIESLKKQ FQASKIELKE EILSKNKVEK DYSNDKRYFL EREQYLLITM ESMESGKLEH EQQTEIMNKD VEELKIRVKE LEATDLSEKN IELERQVGKL TAQLQDFESS RGYEQELILS QLEVQKRKLD KQLIDLQSKL ERNPELVEFD KLHEELQTER MNMPVQVCYV HVSFHKDIE NLIYRKMNNN K
Source	Mammalian cell
Target Names	Y116A8C.11
Protein Names	Recommended name: Uncharacterized protein Y116A8C.11
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.