



Recombinant Human F-box/LRR-repeat protein 17 (FBXL17)

Product Code	CSB-EP866216HU-B
Abbreviation	FBXL17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UF56
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGHLLSKEPR NRPSQKRPRC CSWCRRRRRPL LRLPRRTPAK VPPQPAAPRS RDCFFRGPCM LCFIVHSPGA PAPAGPEEEP PLSPPPRDGA YAAASSSQHL ARRYAALAAE DCAAAARRFL LSSAAAAAAA AASASSPASC CKELGLAAAA AWEQQGRSLF LASLGPVRFL GPPAAVQLFR GPTPSPAELP TPPEMVCKRK GAGVPACTPC KQPRCGGGGC GGGGGGGGGG GPAGGGASPP RPPDAGCCQA PEQPPQPLCP PPSSPTSEGA PTEAGGDAVR AGGTAPLSAQ QQHECGDADC RESPENPCDC HREPPPETPD INQLPPSILL KIFSNLSLDE RCLSASLVCK YWRDLCLDFQ FWKQLDLSSR QQVTDELLEK IASRSQNIIE INISDCRSMS DNGVCVLA FK CPGLLRYTAY RCKQLSDTSI IAVASHCPLL QKVHVGNQDK LTDEGLKQLG SKCRELKDIH FGQCYKISDE GMIVIAKGCL KLQRIYMQEN KLVTDQSVKA FAEHCPELQY VGFMGCSVTS KGVIHLTKLR NLSSLDLRHI TELDNETVME IVKRCKNLSS LNLCLNWIIN DRCVEVIAKE GQNLKELYLV SCKITDYALI AIGRYSMTIE TVDVGWCKEI TDQGATLIAQ SSKSLRYLGL MRCDKVNEVT VEQLVQQYPH ITFSTVLQDC KRTLERAYQM GWTPNMSAAS S
Source	E.coli
Target Names	FBXL17
Protein Names	Recommended name: F-box/LRR-repeat protein 17 Alternative name(s): F-box and leucine-rich repeat protein 17 F-box only protein 13
Expression Region	1-701
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.