



# Recombinant Human F-box/LRR-repeat protein 17 (FBXL17)

<b>Product Code</b>	CSB-MP866216HU
<b>Abbreviation</b>	FBXL17
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UF56
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGHLLSKEPR NRPSQKRPRC CSWCRRRRRPL LRLPRRTPAK VPPQPAAPRS RDCFFRGPCM LCFIVHSPGA PAPAGPEEEP PLSPPPRDGA YAAASSSQHL ARRYAALAAE DCAAAARRFL LSSAAAAAAA AASASSPASC CKELGLAAAA AWEQQGRSLF LASLGPVRF L GPPAAVQLFR GPTPSPAELP TPPEMVCKRK GAGVPACTPC KQPRCGGGGC GGGGGGGGGG GPAGGGASPP RPPDAGCCQA PEQPPQPLCP PPSSPTSEGA PTEAGGDAVR AGGTAPLSAQ QQHECGDADC RESPENPCDC HREPPPETPD INQLPPSILL KIFSNLSLDE RCLSASLVCK YWRDLCLDFQ FWKQLDLSSR QQVTDELLEK IASRSQNIIE INISDCRSMS DNGVCVLA FK CPGLLRYTAY RCKQLSDTSI IAVASHCPLL QKVHVGNQDK LTDEGLKQLG SKCRELKDIH FGQCYKISDE GMIVIAKGCL KLQRIYMQEN KLVTDQSVKA FAEHCPELQY VGFMGCSVTS KGVHILTKLR NLSSLDLRHI TELDNETVME IVKRCKNLSS LNLCLNWIIN DRCVEVIAKE GQNLKELYLV SCKITDYALI AIGRYSMTIE TVDVGWCKEI TDQGATLIAQ SSKSLRYLGL MRCDKVNEVT VEQLVQQYPH ITFSTVLQDC KRTLERAYQM GWTPNMSAAS S
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FBXL17
<b>Protein Names</b>	Recommended name: F-box/LRR-repeat protein 17 Alternative name(s): F-box and leucine-rich repeat protein 17 F-box only protein 13
<b>Expression Region</b>	1-701
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.