



Recombinant Human Inositol polyphosphate 5-phosphatase K (INPP5K)

Product Code	CSB-BP866277HU
Abbreviation	INPP5K
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BT40
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSRKLSGPK GRRLSIHVVT WNVASAAPPL DLSDLLQLNN RNLNLDIYVI GLQELNSGII SLLSDDAFND SWSSFLMDVL SPLSFIKVSH VRMQGILLV FAKYQHLPYI QILSTKSTPT GLFGYWGNKG GVNICKLYG YYVSIINCHL PPHISNNYQR LEHFDRIEM QNCEGRDIPN ILDHDLIWF GDMNFRIEDF GLHFVRESIK NRCYGGLEWEK DQLSIAKKHD PLLREFQEGR LLFPPTYKFD RNSNDYDTSE KKRKPAWTDR ILWRLKRQPC AGPDTPIPPA SHFSLSLRGY SSHMTYGISD HKPVSGETFDL ELKPLVSAPL IVLMPEDLWT VENDMMVSY STSDFPSSPW DWIGLYKVGL RDVNDYVSYA WVGDSKVSCS DNLNQVYIDI SNIPTTEDEF LLCYYSNSLR SVVGISRPFQ IPPGSLREDP LGEAQPQI
Source	Baculovirus
Target Names	INPP5K
Protein Names	Recommended name: Inositol polyphosphate 5-phosphatase K EC= 3.1.3.56 Alternative name(s): Skeletal muscle and kidney-enriched inositol phosphatase
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein with 5-phosphatase activity toward polyphosphate inositol. The protein localizes to the cytosol in regions lacking actin stress fibers. It is thought that this protein may negatively regulate the actin cytoskeleton. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.