



# Recombinant Human THUMP domain-containing protein 2 (THUMPD2)

<b>Product Code</b>	CSB-EP866279HU-B
<b>Abbreviation</b>	THUMPD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BTF0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSEARGEPGS GPEAGARFFC TAGRGLEPFV MREVRARLAA TQVEYISGKV FFTTCSDLNM LKKLKSARL FLLIKKQFPL IISSVSKGKI FNEMQRLINE DPGSQLNAIS IWKNLLELDA KKEKLSQRDD NQLKRKVGGEN EIIAKKLIKIE QMOKIEENRD CQLEKQIKKEE TLEQRDFTTK SEKFQEEEFQ NDIEKAIDTH NQNDLTFRVS CRCSGTIGKA FTAQEVGKVI GIAIMKHFGW KADLRNPQLE IFIHLNDIYS VVGIPVFRVS LASRAYIKTA GLRSTIAWAM ASLADIKAGA FVLDPMCGLG TILLEAAKEW PDVYYVGADV SDSQLLGTWD NLKAAGLEDK IELLKISVIE LPLPSESVDI IISDIPFGKK FKLKDKIKSI LQEMERVLHV GGTIVLLLSE DHHRRLTDCK ESNIPFNSKD SHTDEPGIKK CLNPEEKTGA FKTASTSFEA SNHKFLDRMS PFGSLVPVEC YKVS LGKTD A FICKYKKS HS SGL
<b>Source</b>	E.coli
<b>Target Names</b>	THUMPD2
<b>Protein Names</b>	Recommended name: THUMP domain-containing protein 2
<b>Expression Region</b>	1-503
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.