



Recombinant Human WD repeat-containing protein 85 (WDR85)

Product Code	CSB-BP866281HU
Abbreviation	WDR85
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BTV6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMGCFALQTV DTELTADSVE WCPLQGCRHL LACGTYQLRR PEDRPAGPQN KGGMEVKEPQ VRLGRLFLYS FNDNNSIHPL VEVQRKDTSA ILDMKWCHIP VAGHALLGLA DASGSIQLLR LVESEKSHVL EPLSSLALEE QCLALSLDWS TGKTGRAGDQ PLKIISSDST GQLHLLMVNE TRPRLQKVAS WQAHQFEAWI AAFNYWHPEI VYSGGDDGLL RGWDTRVPGK FLFTSKRHTM GVCSIQSSPH REHILATGSY DEHILLWDTR NMKQPLADTP VQGGVWRIKW HPFHHLHLLA ACMHSGFKIL NCQKAMEERQ EATVLTSHTL PDSLVIYADW SWLLFRSLQR APSWSFPSNL GTKTADLKGA SELPTPCHEC REDNDGEGHA RPQSGMKPLT EGMRKNGTWL QATAATTRDC GVNPEEADSA FLLATCSFY DHALHLWEWE GN
Source	Baculovirus
Target Names	DPH7
Protein Names	Recommended name: WD repeat-containing protein 85
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.