



# Recombinant Human WD repeat-containing protein 85 (WDR85)

<b>Product Code</b>	CSB-MP866281HU
<b>Abbreviation</b>	WDR85
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BTV6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMGCFALQTV DTELTADSVE WCPLQGCRHL LACGTYQLRR PEDRPAGPQN KGGMEVKEPQ VRLGRLFLYS FNDNNSIHPL VEVQRKDTSA ILDMKWCHIP VAGHALLGLA DASGSIQLLR LVESEKSHVL EPLSSLALEE QCLALSLDWS TGKTGRAGDQ PLKIISSDST GQLHLLMVNE TRPRLQKVAS WQAHQFEAWI AAFNYWHPEI VYSGGDDGLL RGWDTRVPGK FLFTSKRHTM GVCSIQSSPH REHILATGSY DEHILLWDTR NMKQPLADTP VQGGVWRIKW HPFHHLHLLA ACMHSGFKIL NCQKAMEERQ EATVLTSHTL PDSLVIYADW SWLLFRSLQR APSWSFPSNL GTKTADLKGA SELPTPCHEC REDNDGEGHA RPQSGMKPLT EGMRKNGTWL QATAATTRDC GVNPEEADSA FLLATCSFY DHALHLWEWE GN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	DPH7
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 85
<b>Expression Region</b>	1-452
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.