



Recombinant Human Coiled-coil domain-containing protein 28B (CCDC28B)

Product Code	CSB-MP866287HU
Abbreviation	CCDC28B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BUN5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDDKKKKRSP KPCLAQPAQA PGTLRRVPVP TSHSGSLALG LPHLPSPKQR AKFKRVGKEK CRPVLGGGS GSAGTPLQHS FLTEVTDVYE MEGGLNLLN DFHSGRLQAF GKECSFEQLE HVREMKEKLA RLHFSLDVCG EEEDDEEEED GVTEGLPEEQ KKTMAADRNLQ QLLSNLEDLS NSIQKLHLAE NAEPEEQSAA
Source	Mammalian cell
Target Names	CCDC28B
Protein Names	Recommended name: Coiled-coil domain-containing protein 28B
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The product of this gene localizes to centrosomes and basal bodies. It interacts and colocalizes with several proteins associated with Bardet-Biedl syndrome (BBS).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.