



Recombinant Human Ferritin heavy polypeptide-like 17 (FTHL17)

Product Code	CSB-MP866310HU
Abbreviation	FTHL17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXU8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MATAQPSQVR QKYDTNCDAA INSHITLELY TSYLYLSMAF YFNRDDVALE NFFRYFLRLS DDKMEHAQKL MRLQNLRGGH ICLHDIRKPE CQGWESGLVA MESAFHLEKN VNQSLLDLYQ LAVEKGDPQL CHFLESHYLH EQVKTIKELG GYVSNLRKIC SPEAGLAEYL FDKLTLGGRV KET
Source	Mammalian cell
Target Names	FTHL17
Protein Names	Recommended name: Ferritin heavy polypeptide-like 17 Alternative name(s): Cancer/testis antigen 38 Short name= CT38
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.