



Recombinant Human Leucine-rich repeat-containing protein 2 (LRRC2)

Product Code	CSB-MP866321HU
Abbreviation	LRRC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYS8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGHKVVVFDI SVIRALWETR VKKHKAQWQKK EVERLEKSAL EKIKEEWNFV AECRRKGIPQ AVYCKNGFID TSVRLLDKIE RNTLTRQSSL PKDRGKRSSA FVFELSGEHW TELPDSLKEQ THLREWYISN TLIQIPTYI QLFQAMRILD LPKNQISHLP AEIGCLKNLK ELNVGFNYLK SIPPELGDCE NLERLDCSGN LELMELPFEL SNLKQVTFVD ISANKFSSVP ICVLRMSNLQ WLDISSNLT DLPQDIDRLE ELQSFLLYKN KLTYPYSML NLKKLTLVVG SGDHLVELPT ALCDSSTPLK FVSLMDNPID NAQCEDGNEI MESERDRQHF DKEVMKAYIE DLKERESVPS YTTKVSFSLQ L
Source	Mammalian cell
Target Names	LRRC2
Protein Names	Recommended name: Leucine-rich repeat-containing protein 2
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.