



Recombinant Human Tripartite motif-containing protein 4 (TRIM4)

Product Code	CSB-MP866336HU
Abbreviation	TRIM4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C037
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEAEDIQEEL TCPICLDYFQ DPVSIIEGHN FCRGCLHRNW APGGGPFPCP ECRHPSAPAA LRPNWALARL TEKTQRRRLG PVPPGLCGRH WEPLRLFCED DQRPVCLVCR ESQEHQTHAM APIDEAFESY RTGNFDIHVD EWKRRLIRLL LYHFKQEEKL LKSQRNLVAK MKKVMHLQDV EVKNATQWKD KIKSQRMRIS TEFSKLHNFL VEEEDLFLQR LNKEEEETKK KLNENTLKLN QTIASLKCLI LEVGEKSQAP TLELLQNPKE VLTRSEIQDV NYSLEAVKVK TVCQIPLMKE MLKRFQVAVN LAEDTAHPKL VFSQEGRYVK NTASASSWPV FSSAWNYFAG WRNPQKTAHV ERFQHLPCVL GKNVFTSGKH YWEVESRDSL EVAVGVCRED VMGITDRSKM SPDVGWAIY WSAAGYWPLI GFPGTPTQQE PALHRVGVYL DRGTGNVSFY SAVDGVHLHT FSCSSVSRLR PFFWLSPLAS LVIPPVTDRK
Source	Mammalian cell
Target Names	TRIM4
Protein Names	Recommended name: Tripartite motif-containing protein 4 Alternative name(s): RING finger protein 87
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies. Its function has not been identified. Alternative splicing of this gene generates two transcript variants.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.