



Recombinant Human BMP/retinoic acid-inducible neural-specific protein 2 (BRINP2)

Product Code	CSB-MP866338HU
Abbreviation	BRINP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0B6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	VSATAAA VVPEQHASVA GQHPLDWLLT DRGPFHRAQE YADFMERYRQ GFTTRYRIYR EFARWKVNNL ALERKDFFSL PLPLAPEFIR NIRLLGRRPN LQQVTENLIK KYGTHFLLSA TLGGEESLTI FVDKQKLGRK TETTGGASII GGSGNSTAVS LETLHQLAAS YFIDRESTLR RLHHIQIATG AIKVTETRTG PLGCSNYDNL DSVSSVLVQS PENKVQLLGL QVLLPEYLRE RFVAAALSII TCSSEGELVC KENDCWCKCS PTFPECNCPD ADIQAMEDSL LQIQDSWATH NRQFESEEF QALLKRLPDD RFLNSTAISQ FWAMDTSLQH RYQQLGAGLK VLFKKTHRIL RRLFNLCRC HRQPRFRLPK ERSLSYWWNR IQSLLYCGES TFPGTFLEQS HSCTCPYDQS SCQGPIPCAL GEGPACAHCA PDNSTRCGSC NPGYVLAQGL CRPEVAESLE NFLGLETDLQ DLELKYLLQK QDSRIEVHSI FISNDMRLGS WFDPSWRKRM LLTLKSNKYK PGLVHVMLAL SLQICLTKNS TLEPVMAIYV NPFGGSHSES WFMPVNEGSF PDWERTNVDA AAQCQNWTIT LGNRWKTTFFE TVHVYLRSLR KSLDDSSNET IYYELEMEDT PSKNLGYMKI NTLQVFGYSL PFDPAIRDL ILQLDYPYTQ GSQDSALLQL IELRDRVNQL SPPGKVRDL FSCLLRHRLK LANNEVGRIQ SSLRAFNSKL PNPVEYETGK LCS
Source	Mammalian cell
Target Names	BRINP2
Protein Names	Recommended name: Protein FAM5B Alternative name(s): BMP/retinoic acid-inducible neural-specific protein 2 DBCCR1-like protein 2
Expression Region	34-783
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.