



Recombinant Arabidopsis thaliana Eukaryotic translation initiation factor 3 subunit H (TIF3H1)

Product Code	CSB-MP866372DOA
Abbreviation	TIF3H1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C5Z2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATMARSFLQ AISKDEAVAP PLRVVQIEGL AVLKIIKHCK EFSPTLVTGQ LLGLDVGSVL EVTNCFFFPV RDDDEEIEAD GANYQLEMMR CLREVNVDNN TVGWYQSTVL GSYQTVELIE TFMNYQENIK RCVCIIYDPS KADLGVLALK ALKLSDSFME LYRGGNFTGE KLREKNFSWM DIFEEIPIKV SNSALVSAFM TELETDTPVS QGDYDRLHSS TTPFLENME FLIKCDDLS MEQQKFQYYY RNLSRQQAQQ QAWLQKRRTN NMARKSAGEE PLPEEDPSNP IFKAIPEPSR LESFLITNQV SNFCGQINGV AGQNFSRLYL TKALHDN
Source	Mammalian cell
Target Names	TIF3H1
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit H Short name= eIF3h Alternative name(s): Eukaryotic translation initiation factor 3 subunit 3 eIF-3-gamma eIF3 p38 subunit
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.