



Recombinant Arabidopsis thaliana Gibberellin 2-beta-dioxygenase 7 (GA2OX7)

Product Code	CSB-EP866377DOA-B
Abbreviation	GA2OX7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C6I4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASQPPFKTN FCSIFGSSFP NSTSESNTNT STIQTSGIKL PVIDLSHLTS GEEVKRKRCV KQMVA-AAKEW GFFQIVNHGI PKDVFEMMLL EEKCLFDQPF SVKVRERFSD LSKNSYRWGN PSATSPAQYS VSEAFHIILS EVSRISDDRN NLRTIVETYV QEIARVAQMI CEILGKQVNV SSEYFENIFE LENSFLRLNK YHPSVFGSEV FGLVPHTDTS FLTILSQDQI GGLELENNGQ WISVKPCLEA LTVNIGDMFQ ALSNGVYQSV RHRVISPANI ERMSIAFFVC PYLETEIDCF GYPKKYRRFS FREYKEQSEH DVKETGDKVG LSRFLI
Source	E.coli
Target Names	GA2OX7
Protein Names	Recommended name: Gibberellin 2-beta-dioxygenase 7 EC= 1.14.11.13 Alternative name(s): GA 2-oxidase 7 Gibberellin 2-beta-hydroxylase 7 Gibberellin 2-oxidase 7
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.