



Recombinant Arabidopsis thaliana Probable F-box protein At1g44080 (At1g44080)

Product Code	CSB-MP866381DOA
Abbreviation	At1g44080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C6X9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGKVGWSDLH EDLIDLLANN LSSNINLLRF RSICKPWRST VATKKRLHNNH FERNLPTFFK KKTVVSPSTF FRVTLPSPCR NKGWLIKNRQ VSESSKNLL SPLSGKTITP SDKTLDLLKV ECFRDSSILQ LFADSDRVVF LDNVFFVDF KNEIWCKSG EETRHWRTRIN NEEAKGFLDI ILHKGKIYAL DLTGAIWWIS LSELSIQYG PSTPVDFYEI DNCCKEKLVE YCGELCVVHR FYKKFCVKRV LTERTVCFKV YKMDKNLVEW VEVSSLGDKA LIVATDNCFL VLASEYYGCL ENAIYFNDRE DVSVFKLDDC SIITKQINSS SSSSHQSCFQ MFSTSF
Source	Mammalian cell
Target Names	At1g44080
Protein Names	Recommended name: Probable F-box protein At1g44080
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.