



# Recombinant Arabidopsis thaliana Cytosolic sulfotransferase 18 (SOT18)

<b>Product Code</b>	CSB-BP866401DOA
<b>Abbreviation</b>	SOT18
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C9C9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESETLTAKA TITTTTLPSH DETKTESTEF EKNQKRYQDL ISTFPHEKGW RPKEPLIEYG GYWWLPSLLE GCIHAQEFFFQ ARPSDFLVCS YPKTGTTWLK ALTFAIANRS RFDDSSNPLL KRNPFHFVPY IEIDFPFFPE VDLVKDKGNT LFSTHIPYEL LPDSVVKSGC KMVYIWREPK DTFISMWTFH HKERTELGPV SNLEESFDMF CRGLSGYGPY LNHILAYWKA YQENPDRILF LKYETMRADP LPYVKSLAEF MGHGFTAEFE EKGVVVEKVVN LCSFETLKNL EANKGEKDRE DRPGVYANSA YFRKGVGDW SNYLTPEMAA RIDGLMEEKF KGTGLLEHGK
<b>Source</b>	Baculovirus
<b>Target Names</b>	SOT18
<b>Protein Names</b>	Recommended name: Cytosolic sulfotransferase 18 Short name= AtSOT18 EC= 2.8.2.- Alternative name(s): Desulfo-glucosinolate sulfotransferase B Sulfotransferase 5B Short name= AtST5b
<b>Expression Region</b>	1-350
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.