



Recombinant Arabidopsis thaliana Nudix hydrolase 26, chloroplastic (NUDT26)

| | |
|--------------------------|---|
| Product Code | CSB-YP866413DOA |
| Abbreviation | NUDT26 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9CAF2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | SSMESPP EGYRRNVGVC LMNSSKKIFT ASRLDIPSAW QMPQGGIDEG EDPRVAVMRE LKEETGVHSA EILAEAPHWI TYDFPPDVRE KLKVRWGSDW KGQAQKWFL KFTGKDEEIN LLGDGTEKPE FGEWSWTSPD QVVENAVEFK KPVYKEVMSA FASHLQ |
| Source | Yeast |
| Target Names | NUDT26 |
| Protein Names | Recommended name: Nudix hydrolase 26, chloroplastic Short name= AtNUDT26 EC= 3.6.1.- Alternative name(s): Bis(5'-nucleosyl)-tetrphosphatase (asymmetrical) EC= 3.6.1.17 |
| Expression Region | 54-216 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |