



Recombinant Arabidopsis thaliana Ras-related protein RABF1 (RABF1)

Product Code	CSB-BP866419DOA
Abbreviation	RABF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CB01
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	GCASSLPDR NSGTLGSLSN SENAVPADAK NLRVKLVLLG DSGVKGKSCIV LRFVVRGQFDA TSKVTVGASF LSQTIALQDS TTVKFEIWDT AGQERYSALA PLYYRGAGVA VIVYDITSPE SFKKAQYWVK ELQKHGSPDI VMALVGNKAD LHEKREVPTE DGMELAEKNG MFFIETSAKT ADNINQLFEE IGRKLRPAP SS
Source	Baculovirus
Target Names	RABF1
Protein Names	Recommended name: Ras-related protein RABF1 Short name= AtRABF1 Alternative name(s): Ras-related protein Ara-6 Ras-related protein Rab5C Short name= AtRab5C
Expression Region	2-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.