



# Recombinant *Lactococcus lactis* subsp. *lactis* Uncharacterized RNA methyltransferase yljE (yljE)

<b>Product Code</b>	CSB-EP866471LNG-B
<b>Abbreviation</b>	yljE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9CGB9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Lactococcus lactis</i> subsp. <i>lactis</i> (strain IL1403) ( <i>Streptococcus lactis</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCYNMDMTNF KKNDIFEAEV LDLTHEGQGV VKIDSFPFFV DNALPGERIK MHVLKVGKSF GFGRVDEFIS QSKHRVTQLN LDYLRGTGIAD FGHLSYDEQV KFKHKQVVDI LRKTAGKSDY PVLPPVEAIK TINYRNKAQI PVQMVDGLLT TGFYRKNSTH LIPVEDFYIQ HPEIDAVVLF LRDEFKINL QAYDEKTRKG WLRNIVVRRR FHTGEMMITL VVTSKKLPEK VELVIDHLVE RFTNIKSVQL NINSGTGSFI LGKEFILLYG KDFITDMLD KTYQISAPAF YQVNTPPAEK LYETAYEFAG LKPGDVVIDA YSGIGTIGIS MADRVAKVYG MEVVPAAVEN AKRNAQLNEL ENTHYEVGTA EKIMPKWLSE GIKPDVIFVD PPRKGLDEGF IKAAAMMNPR RIVYISCNPA TFARDVVRFE ADGYLLDKVQ PVDLFPQTHH IELVASFNKK
<b>Source</b>	E.coli
<b>Target Names</b>	yljE
<b>Protein Names</b>	Recommended name: Uncharacterized RNA methyltransferase yljE EC= 2.1.1.-
<b>Expression Region</b>	1-460
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.