



# Recombinant Mouse NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial (Ndufs3)

<b>Product Code</b>	CSB-BP866759MO
<b>Abbreviation</b>	Ndufs3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9DCT2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ESAAA DKRPTVRPRS DVTHKQLSAF GEYVAEILPK YVQQVQV SCL DELEICHPD GVIPTLTFLR DHTNAQFKSL ADLTAVDVPT RQNRFEIVYN LLSLRFNSRI RVKTYADELT PIDSIVSVHI AANWYEREVW DMFGVFFFNH PDLRRILTDY GFEGHPFRKD FPLTGYVELR YDDEVKRVVA EPVELAQEFR KFDLNSPWEA FPAYRQPPES LKLEAGDKKP ETK
<b>Source</b>	Baculovirus
<b>Target Names</b>	Ndufs3
<b>Protein Names</b>	Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-30kD Short name= CI-30kD NADH-ubiquinone oxidoreductase 30 kDa subunit
<b>Expression Region</b>	36-263
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.