



Recombinant Mouse NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial (Ndufs3)

Product Code	CSB-EP866759MO-B
Abbreviation	Ndufs3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DCT2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	ESAAA DKRPTVRPRS DVTHKQLSAF GEYVAEILPK YVQQVQVSCL DELEICHPD GVIPTLTFLR DHTNAQFKSL ADLTAVDVPT RQNRFEIVYN LLSLRFNSRI RVKTYADELT PIDSIVSVHI AANWYEREVW DMFGVFFFNH PDLRRILTDY GFEGHPFRKD FPLTGYVELR YDDEVKRVVA EPVELAQEFR KFDLNSPWEA FPAYRQPPES LKLEAGDKKP ETK
Source	E.coli
Target Names	Ndufs3
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-30kD Short name= CI-30kD NADH-ubiquinone oxidoreductase 30 kDa subunit
Expression Region	36-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.