



Recombinant Danio rerio Transcription factor Sox-19b (sox19b)

Product Code	CSB-MP866763DIL
Abbreviation	sox19b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DDD7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MYSMMEHELK TAGPPHTLQH SPGMSPPGSG VGNAHHVSKT ACPPGVDPMD KVKRPMNAFM VWSRGQRRKM AQENPKMHNS EISKRLGAEW KLLTDVEKRP FIDEAKRLRA VHMKEYPDYK YKPRRKTKAL MKKDNPVGKY PLAAGNLLAS AVAQGGGSP RMDGYGWGHA GGYMGMQGDA LGYPQQLHRY DLSALQYPAA QPYMNSASSY SQMSYSSSPQ QPSPVMSMVK PEPLSHSPTG VPNHHRGAFQ GDLRDMISMY IPGGDTSESS NQRAYPGVQQ HYLGGTVPLT HI
Source	Mammalian cell
Target Names	sox19b
Protein Names	Recommended name: Transcription factor Sox-19b Alternative name(s): Protein sox-31
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.