



Recombinant *Oryzias javanicus* G2/mitotic-specific cyclin-B2 (ccnb2)

Product Code	CSB-BP866775OCAE
Abbreviation	ccnb2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DG99
Product Type	Recombinant Protein
Immunogen Species	<i>Oryzias javanicus</i> (Javanese ricefish) (<i>Aplocheilus javanicus</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSSVEAVTQQ QLLAAENPVR MGKGAADPRR AALGEITNRN AAAAANRKMG PSKKQPKPSC AQKPQPVVHT SAGDPAPISA DMSMKVEQDL SQAFSEVLML AVQDVDEQDA DQPQLCSQYV KDIYKYLHTL EEQQAIRPNY MQGYSVTEHM RALLVDWL VQ VHSRFQLLQE TLYLTVAILD RFLQVHPVSR RKLQLVGVTA MLVACKYEEM YPPEVGDFAY ITDDAFTKFQ IVEMEQVILR SLGFQLGRPL PLHFLRRASK VADADVEKHT LAKYLLELTL LDYHMHVHYRP SEAAAAALCL SQLLLDGLPW SLEQQHYSTY DEQHLKPI MQ LMAKNVVQVT EGRTKFLAVK KKYSSSKLMK ISLIPQLNSS TIKVMAEALQ NP
Source	Baculovirus
Target Names	ccnb2
Protein Names	Recommended name: G2/mitotic-specific cyclin-B2
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.