



Recombinant Mouse Peroxisomal acyl-coenzyme A oxidase 3 (Acox3)

Product Code	CSB-MP866807MO
Abbreviation	Acox3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EPL9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	GSLPEEKDS ALWSDTPKGP LSAYRARASF NSGELLLFWD GQDVIHFKKT IFSTLENDPL FARSYGADLP LEKLRELNFL RCKRVFEYGF FKVEELLKNP LKILVLINCL GMYDWSLANK CVLHMLVFGT TVFVSGSEKH FKYLEKIYSL EIFGCFALTE LSHGSNTKAM RTTAHYDPDT QEFILHSPDF EAAKFWVGNL GKTATHAVVF AQLYMPDGQC HGLHSFLVQI RDTKTLLPMT GVMVGDIGKK LGQNGLDNGF AMFNKVRIPR QNLLDRTGNI TSEGTYNSPF KDVRQRLGAS LGSLSSGRIS IISMSVVNLK LAVSIAIRFS ATRCQFGPTD KEEIPVLEYP LQQWRILPYL AAAYALDHFS KTIFMDLIEV QSARLRGDHS DQQAELGREI HALASAGKPL ASWTAQRGIQ ECREACGGHG YLAMNRFGLD RNDNDPNCTY EGDNNVLLQQ TSNYLLSLE PPLQDGAHFT SPLKTVDFLE AYPGILGQKF LGSSKADWMD SAAPLAAYRW LVCYLLQESH RRYCQEKKSR GSDFEARNNS QVYGCRPLAL AFMELTVMQR FHEHIHSSGL SPSLRTVLGR LSTLYGLWCL SQHMALLYRG GYISGEQTGR AMEDAILTLC EQLKDDAVAL VDVIAPSDFV LNSPIAKADG ELYKNLWAAV LQQNGVLERA AWWPEFSANK SVADRLKSQL
Source	Mammalian cell
Target Names	Acox3
Protein Names	Recommended name: Peroxisomal acyl-coenzyme A oxidase 3 EC= 1.3.3.6 Alternative name(s): Branched-chain acyl-CoA oxidase Short name= BRCACox Pristanoyl-CoA oxidase
Expression Region	2-700
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Acyl-Coenzyme A oxidase 3 also know as pristanoyl -CoA oxidase (ACOX3)is involved in the desaturation of 2-methyl branched fatty acids in peroxisomes.



Unlike the rat homolog, the human gene is expressed in very low amounts in liver such that its mRNA was undetectable by routine Northern-blot analysis or its product by immunoblotting or by enzyme activity measurements. However the human cDNA encoding a 700 amino acid protein with a peroxisomal targeting C-terminal tripeptide S-K-L was isolated and is thought to be expressed under special conditions such as specific developmental stages or in a tissue specific manner in tissues that have not yet been examined.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.