



Recombinant Mouse Fez family zinc finger protein 2 (Fezf2)

Product Code	CSB-BP866847MO
Abbreviation	Fezf2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ESP5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MASSASLETM VPPACPRAGA SPATSKTLAF SIERIMAKTS EPRAPFEPRP AALEADSSQS KKLLNLCSP L PCMIPLQPLG YEVPSTLLS YSEFWKSSLR AGGGGGGGSG GGAPVCGASG LCKTNCGVCC KAEGLAPSA LPAGRVIKPQ VINQAVGLPA SGSLYYFNYL DSTAYPPSEL LGGHLFPSGL LNAQAPTSLA AHPKLFLEN AKLASLAADK FPHPASYPHK ERLHAPLEQV LKENSALTAE RGGVKSHSKL PGGSTDSKPK NFTCEVCGKV FNAHYNLTRH MPVHTGARP VCKVCGKGFR QASTLCRHKI IHTQEKPHKC NQCGKAFNRS STLNTHIRIH AGYKPFVCEF CGKGFHQKGN YKNHKLTHSG EKQYKCTICN KAFHQVYNLT FHMHTHNDKK PFTCATCGKG FCRNFDLKKH VRKLHDSVGP TATPSAKDLA RTVQS
Source	Baculovirus
Target Names	Fezf2
Protein Names	Recommended name: Fez family zinc finger protein 2 Alternative name(s): Forebrain embryonic zinc finger-like protein 2 Zinc finger protein 312 Zinc finger protein Fez-like
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.