



Recombinant Arabidopsis thaliana Ubiquitin-conjugating enzyme E2 22 (UBC22)

Product Code	CSB-MP866905DOA
Abbreviation	UBC22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FF66
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MASNENLPPN VIKQLAKELK SLDESPPDGI KVVVNDEDFS QICADIEGPV GTPYENGLFR MKLALSHDFP HSPPKG YFMT KIFHPNVASN GEICVNTLKK DWNPSLGLRH VLSVVRCLLI EPFPESALNE QAGKM LLENY DEYARHARLY TGIHAKPKPK FKTGAISEST TALNVGQTNN ETPGAATAIP SSMTDIKRVT TSAQDQQHVA NVVVAASASV VTTTQKREAG LAKVQADKKK VDARKKSLKR L
Source	Mammalian cell
Target Names	UBC22
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 22 EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein 22
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.