



Recombinant Arabidopsis thaliana Putative B3 domain-containing protein At5g66980 (At5g66980)

Product Code	CSB-EP866918DOA-B
Abbreviation	At5g66980
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGD2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVDLSKDLQF FKVFLPEFGS HELVIPPAPI DMLEKPLPKE AFLVDEIGRL WCVETKTEDT EERFCVFFTK GWQSFANDQS LEFGDFLVFS YDGDSRFSVT IFANDGCKKD VGVVSTTDRS RVSLDEEEDPD DFTKPDRMR DCDCGQSINR KRKRDSVNED PHLVDDKPE YVSTYKTKPE HSEKTQRTVN RAGDTCDISW FPEKKHNGFE ESVYKPKHPH FVRNITRGSL QKLELPLTFL RSNGIELEED IELCDESGKK WPLKILNHDR GFKFSHESWL CFCKSHEMIL TNKCLFEFIV PSNGRCSEIL VRIVSGRLPT TMTKNNYQVI DMQP
Source	E.coli
Target Names	At5g66980
Protein Names	Recommended name: Putative B3 domain-containing protein At5g66980
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.